

WWW.CHISON.COM



ECO 6

More Advanced **More** Powerful

CHISON
Value Beyond Imaging

Easy to operate Easy to go



ONE KEY to selected the desired application with automatic preset



ONE KEY to optimize the image automatically



ONE KEY to save still image and CINE loop



ONE KEY to measure



ONE KEY to make report





ONE KEY to print to PC or video printer (see page 4 option)



More advanced

Continuous Wave (CW) Doppler

Detect blood flow with high velocity and help doctors diagnose with more clinical information.




Tissue Doppler Imaging (TDI)

Tissue Doppler imaging is a novel echocardiography technique that directly measures myocardial velocity. Systolic TD measurements assess left and right ventricular myocardial contractile function. Diastolic TD values reflect myocardial relaxation.

Triplex

Triplex combines B, color, with PW to help the user make diagnosis more conveniently and more accurately, help to improve the efficiency



Auto IMT Function

Automatically traces the intima, and measures the thickness of the intima. This allows you to measure the intima faster, more easily and more accurately.



Without Super Needle
Needle is Not clear



With Super Needle
Needle is Better Visualized

Super Needle

With Super Needle, clinicians can see needle inside tissue more clearly during medical procedures. Needle angle up to $\pm 30^\circ$

SonoBladder (option)

- Automatically measure bladder volume
- Clearly visible in dual format



Auto Thyroid (option)

- Automatically detect thyroid lesion
- Provide the size
- Efficient for the diagnose

AutoBreast (option)


- Automatically detect breast lesion
- Provide the size
- Efficient for the diagnose




More Powerful

- ✓ Abdomen, OB/GYN, Urology, cardiac
- ✓ Phased array probe
- ✓ CW, TDI, B/BC, super needle
- ✓ 2D steer, triplex, Auto IMT, Automatic pw trace, CPA, DPD






Thyroid, CPA Mode



Carotid, 2D steer



Carotid, Triplex



Liver, B Mode




Fetal Head, B Mode




Kidney, B/BC Mode



Cardiac, B/BC Mode



Liver, PW Mode



Cardiac, PW Mode



Convex
C3-A



Linear
L7M-A



Linear
L7S-A



Transvaginal
V5-A



Trans-Rectal
RT-A



Phased Array
P3-A



Micro-Convex
MC3-A



Micro-Convex
MC6-A

CHISON Medical Technologies Co., Ltd.

Sales & Service Contact Address: No.3, Changjiang South Road, Xinwu District, Wuxi, Jiangsu, China 214028

TEL : 0086-510-85310937 **FAX :** 0086-510-85310726 **EMAIL :** export@chison.com.cn

We reserve the right to make changes to this catalogue without prior notice. All rights reserved.
Please contact our local dealer for the latest information.